An Introduction to LabVIEW for 4th year projects



Stephan Eisenhardt, University of Edinburgh

S.Eisenhardt@ed.ac.uk






Objectives

- □ To acquire familiarity with the LabVIEW Programming language
- To be able to write LabVIEW programmes incorporating pre-written and new code
- □ To be able to understand, adopt and modify third-party code
- **To learn how to investigate about unknown functionality**
- □ To acquire the skills needed to complete the 4th year project

Stage 0 : Concepts & Principles

The LabVIEW Concept I


- LabVIEW (in short 'LV') is a graphical programming language developed and marketed by National Instruments:
 - The look and feel is very different from textual programming languages



- All the usual programming concepts are available
- But what is the benefit of reinventing the wheel, causing overhead, another learning curve and paying substantial licence fees?

The LabVIEW Concept II

- LabVIEW is designed to build graphical user interfaces (GUIs) for laboratory instrumentation – also called 'Virtual Instruments' (VIs):
 - A VIs can read out to a single sensor, emulate the front-face of commercial devices or control large integrated systems
 - The idea is to provide a user-friendly interface, tailored to the needs of the application, to an otherwise possibly obscure piece of hardware



- Like with a dashboard of a car, the user interface is what the operator experiences for steering and monitoring of the hardware, while the details of the implementation and the interfaces to the devices are hidden under the bonnet
- Programming LabVIEW is like building a car, running a VI like driving one

Introduction to LabVIEW

Stephan Eisenhardt

Principle of 'Data Flow'

- □ Programme execution follows the principle of 'data flow':
 - Each instruction comes with an interface of input and output parameters, called 'terminals' (e.g. the '+'-operation has 'x' and 'y' as input and the 'sum' as output terminal)
 - An instruction is only executed once *all* its input terminals hold valid data (i.e. the instructions it depends on, all have produced valid data)
 - If there are several instructions which could be executed in parallel, LV has reasonable defaults to chose; but if the order matters, e.g. due to 'race conditions', there are ways to control this
 - Each VI has it own set of terminals and can be called as an instruction in another VI (e.g. also function calls provided by LV are VIs themselves)

Implemented Design Principles

- Encapsulation: (main design feature)
 - Each VI has its user-definable interface, and can be operated as 'black box' element within other VIs once it is working reliably
- □ Type checking: (main design feature)
 - Checked at time of coding, VI cannot run if output and input between two nodes do not match in type
- □ Templates: (now fully integrated since v8.0 and promoted as default)
 - Called 'Express VIs': configurable, pre-written code for standard tasks
- □ Scope of parameters: (supported)
 - VIs run with their own set of local variables, and one can manage different sets of global variables...
- Recursive functions: (supported)
 - VIs can be configured to run as parallel instances and can call themselves
- □ Inheritance: (not in the way you know it from C++, only while coding)
 - Configurations of existing elements may be 'inherited' by new elements you place
- □ Function overloading/Polymorphism: (implemented for library functions)

Provided for many standard functions in the library, but hard to code by one-self
 Introduction to LabVIEW
 Stephan Eisenhardt

Stage 1 : The Way around LabVIEW

How to use this Guide

- **Layout**:
 - Proposed instructions are green
 - Items to note are orange
- □ While reading the following pages:
 - Run LabVIEW in parallel and try all the discussed actions for yourself
 - Feel free to branch out from the guide to browse some of the many other possibilities and return to the text at your own pace
 - At this stage you are deliberately *not* given any example.vi the first stage is about learning to find your way around and ways to help yourself
 - A second stage will deal with the functionality you will need to employ during your project

Introduction to LabVIEW

Getting Started

- □ Start LabVIEW: (used here: version 9.0.1)
 - Windows: double-click LV icon or Start menu \rightarrow All Programs
 - Linux/SunOS: type 'labview &'
- □ Startup screen: select 'New: Blank VI'
 - We start from scratch and keep it simple

 A simple data member and a sum will give the opportunity to learn much about the way around LV




 That should give you the tools and techniques to better find the needed information later on


Basics

□ A LabVIEW programme comprises of two types of screens:

 The 'front panel' (FP): this acts as the user interface, with controls and displays
 associated with it is a 'Controls'
 window (Cnt) to drag & drop GUI
 interface elements



- The 'Block Diagram' (BD): (invoke with Ctrl-e if needed) this contains the computational code which handles the data associated with it is a 'Functions' Functions 🔍 Search 🛛 👫 View 🔻 Express window (Ftn) to drag & drop code 7 **9**, hu Input Signal Analysis Output elements ⊳Σ⁺ ∫⊳ ∎≓≣ 08 Sig Manip Exec Control Arith & Compa
 - it is stored in a file <name>.vi



Introduction to LabVIEW

Stephan Eisenhardt

Favorites

User Libraries
 Select a VI...

Mouse Tips

- Basic manipulation techniques:
 - hovering over items selects the focus
 - left-click selects an item
 - double left-click selects to edit an item
 - right-click opens a menu for configuration and digging further

Learn from a Simple Example

- □ A simple example will show the basic functionality: a data member
 - In the 'Controls' window: under the 'Express' tab open the 'Numerical Controls' sub-tab and select the plain 'Numeric Control' GUI

🔍 Search 🛛 👫 View 🕶

User Controls Select a Control.. abc

Text Ctrls

0

LEDs

Num Ind

- drag & drop it to the FP
- notice what happened:
- 1. on the FP you got an interface

with a data display, increment handle and name label

 \rightarrow play around with it: edit the value, the label, use the increment

 \rightarrow notice the tick button whenever a value is still transient

- 2. on the BD you got an input 'terminal' with the same label
 - \rightarrow it represents the memory allocation and displays the data type

Introduction to LabVIEW

Numeric

1.23

ontrols

Express
 L Numeric Controls

Fill Slide

User Controls

Select a Control...

🔍 Search 🛛 👫 View 🔻

0 5 10

Fill Slide

10-5-0-

Pointer Slide

Color Box

Pointer Slide

à

Knoh

Numeric

010

close


open

Element Manipulation I

- □ Lets dig a bit deeper into what can be done to the data members:
 - right-click the input terminal and select 'Representation':

you see the numerical data types LV can handle by rows, in the different precisions:


- floating point
- signed integer
- unsigned integer
- complex



- by default LV chooses 'double precision' = 32-bit floating point
 - \rightarrow note how colour and icon of the terminal changes for the various types
 - → try it out, e.g. to change to 'complex extended precision' and see how the data display on the FP changes

Element Manipulation II

- And on the front panel:
 - right-click the data display and select 'Visible Items': you see now also the other visual items available:
 - \rightarrow try out the radix, select the 'U8' data type and value '10':
 - \rightarrow click the radix and change to 'Hex', 'Octal' and 'Binary'
 - \rightarrow note that you have to chose an integer data type for these settings to be available
 - \rightarrow try to activate the 'Unit Label', you cannot...
 - \rightarrow you have to chose a floating data type first to use a 'Unit Label', use e.g. 'V' as label
 - \rightarrow now you cannot change the radix to 'Hex', 'Octal' or 'Binary' anymore, nor can you change to integer data types...
 - \rightarrow you must clear the 'Unit Label' again by editing (double-click), to enable integer types once more
- that is all due to the data type checking already performed at the coding stage Introduction to LabVIEW



Element Manipulation III

- Dizzy already? There still is more worth to know at this stage:
 - for the terminal or data display find the 'Data Operations' menu and 'Reinitialize to the Default Value': the data display should show '0' again
 - → enter a value<>0 and 'Make Current Value Default'
 - \Rightarrow provided you save now the .vi, next time you call it you will start with the value you have chosen
 - This works fine: but for transparency I strongly recommend to explicitly initialise all terminals to sensible defaults at the start of running, it makes life (debugging) so much easier... – we will come to that...

Element Manipulation IV

- □ And one more on 'states' of elements:
 - for the data display find the 'Advanced' menu and use 'Hide Control'
 - ... oops... you want it back? \rightarrow terminal: 'Show Control', easy
 - another game on the 'state' of our element:
 data display → Advanced → Enabled State → select Disabled & Grayed
 ... hmm... you still can do everything you have done so far...
 yes, because you are in 'editing mode', 'disabling' means to lock the access at run-time. i.e. for the operator
 - Now: all operations on 'states' of elements, e.g. data operations, colouring, animation, formatting, visibility of components, etc., are available at run-time as well; so called 'Property Nodes' allow for their test and manipulation during program execution... imagine ...
 - you find a comprehensive collection of states of an element in the rightclick menu under 'Properties' → go and explore...

Introduction to LabVIEW

Quo Vadis?

□ You are worried because you 'haven't done anything' yet?


- So wrong! You have achieved so much more than in the usual approach of putting 'some lines of code together and run it'...
- You have gasped a taste of the scope of LabVIEW!
- You have gained a set of tools to query and understand the elements of your programme!
- we continue, and you should expect much more from LabVIEW: as each new principle opens a new dimension in the spectrum of possibilities
- don't shy away from the mightiness of LV, curiously explore
- BTW: Vee is the competitor language developed by Agilent, one of the main competitors of National Instruments
 - when you learn LV, you quickly can pick up Vee as well as they are similar
 - but the languages are not compatible...: imagine how the C-compilers of Microsoft and Borland would look like if there wouldn't be the ANSI consortium to standardise the C-language...


Introduction to LabVIEW

Stephan Eisenhardt

Help Yourself

- Here is what you can do to get more information:
 - the 'Context Help' window is extremely useful:
 - · it gives you the synopsis of the element you mouse-over
 - describes all I/O ports of an element
 - gives a link to the full reference (for all library functions)
 - to invoke do either:
 - $Help \rightarrow Show Context Help$
 - push the '?' button
 - Ctrl-h





- for free textual search go to LabVIEW help:
 - Help \rightarrow Search the LabVIEW Help...
 - Ctrl-?



Introduction to LabVIEW

Reading Third-Party Code


- □ Learning from examples is a great way to progress:
 - menu: Help \rightarrow Find examples...
 - explore the possibilities to:
 - narrow search to specific hardware
 - maintain your favourites
 - online information in the LabVIEW Zone
 - read existing third-party code
 - some examples are simple, some complex
 - to make the most of it: use the following strategies:
 - follow the flow of data

understand the I/O of elements and the used functionality, get ideas

• follow the variables

use right-click \rightarrow Find xxx tools to see where data is used/re-used in the code

 use the debugging tools to study the data flow and sequence of events



Introduction to LabVIEW

Stephan Eisenhardt

Ways to Build: the second element

- □ Now relax as this is easy but there are several ways to Rome:
 - revert back to the default settings of the first Numeric Control and add a second Numeric Control
 - which way have you done it?
 - via the Controls window again? OK
 - you also could have done via Copy & Paste (Ctrl-c & Ctrl-v) on either


ļ	the FB	or	the BD \rightarrow	try both, and delete (Ctrl-x) or undo (Ctrl-z)
	Numeric 0		Numeric (123)	the advantage and risk with Copy&Paste is
	Numeric 2		Numeric 2	that you also copy all configurations you have done to the source (did you really reset
ļ				everything?)

- note the different behaviour if you perform the actions on the FP or BD:
 - if you place the new element on the FP: the new element on the BD will be placed in (approximately) the same relative position to the first element
 - if you place the new element on the BD: the new element on the FP will be placed as near as possible (without overlap) to the top-left corner of the FP

Introduction to LabVIEW

Ways to build: chose an operation

- □ As always in LV, again you have several ways to do it:
 - in the Functions window: activate the extension button follow: → Programming → Numeric drag the 'Add' operation to the BD
 - note the tree structure of available functionality



Functions

Measurement I/C
 Instrument I/O
 Vision and Motion
 Mathematics

- go & find in the 'Programming' branch: While Loop, Local Variable, Array Subset, Transpose Matrix, Bundle Cluster, π constant, True constant, type conversion Boolean to (0,1), Carriage Return constant, Max&Min test, Wait Until Next ms Multiple, Three Button Dialog, Scan From File, Property Node, Acquire Semaphore, Play Sound File
- note that many of these occur redundantly again in other parts of the tree
- others do not: so go & find: Gauss Peak Fit, Numeric Integration, Sine Wave Generator, Random Noise Generator, Convolution, Normalisation, 1/f Filter, Fast Fourier Transformation, Initialize Mouse

Ways to build: use the context menu

□ No, that was only one way so far, here comes the second:

Visible Items Find Control Hide Control Change to Indicator Change to Constant Description and Tim... right-click on one input terminal: depending on its data type you will get a palette entry which provides the most common elements this terminal would connect to, and links to further related sub-palettes

Numeric Palette			A	dd		
Create Data Operations	\triangleright	Þ	\blacktriangleright	\Rightarrow	÷R] <u>]32</u>) ▶]DBL
Advanced 🕨 🕨	Add	Subtract	Multiply	Divide	Quotient & R	Conversion
✓ View As Icon	+	₽	Σ	Þ	∃ }	
Representation	Increment	Decrement	Add Array Ele	Multiply Array	Compound Ar	Data Manipul
Properties	Þ			\triangleright	×2	8+iY
	Absolute Value	Round To Ne	Round Towar	Round Towar	Scale By Pow	Complex
	\triangleright	h	D			mx+b ~+^√
	Square Root	Square	Negate	Reciprocal	Sign	Scaling
	123	• Enum	Ring	r de la companya de l	EXPR.	FXP D.CC
	Numeric Cons	Enum Constant	Ring Constant	Random Num	Expression N	Fixed-Point
	+00	-00	E			₩
	+Inf	-Inf	Machine Epsilon			Math Constants


this context-menu provides a tailored selection of options fitting to the element it is used for, in general it is much faster to use this one – and it seems pretty comprehensive

- a look ahead: in the same way under the 'Create' entry you find an exhaustive list of related elements to create
- the clou with this is: the new element is already configured and linked to the 'parent' element, it inherits its properties
- try to create a 'Property Node' to programmatically control the foreground colour of the label...

Ways to build: use inheritance

□ This way of coding, by using context menus, is powerful and efficient:

– now verify the inherited properties:





 right-click the coloured part of the property node → Link To → Pane:

see to which terminal it is linked


- note that the property Label.FGColor is already configured
- right-click the property itself: note the additional entries in the menu, enabling manipulation, e.g. to add an element and Change to Write
- when you left-click on the property you can chose a different one
- now you can get rid of the Property Node again


Introduction to LabVIEW

Wiring Basics

- □ Now we connect the elements:
 - remember that 'Data Flow' regulates the flow of the programme
 - the data travels between elements through 'wires'
 - your Block Diagram should look like this again:
 - mouse-over the terminal and '+' operation:
 - note the changes to the mouse and the icons: the mouse looks like a little spindle when you hover over the I/O ports of the elements – which are indicated by orange circles and black background for the active one
 - by pointing & dragging: string a wire from one input terminal to the xinput of the '+' operation
 - and add the wire for the second
 - note: output goes to the right input comes from the left


Introduction to LabVIEW






More Wire Skills

- You will use this a lot...:
 - mouse over a wire: actually slightly next to it, so that the mouse pointer is a spindle
 - click and start to draw an extension to the wire
 - click again to create fixed support points while you drag
 - right-click to abandon your wiring
 - try again and now double-click to end the wire mid-air
 - note that the new wire, by dashing, is indicated as non-functional
 - also the original wire bears a cross to say that it doesn't work anymore
 - mouse-over the cross and see what is wrong: fair enough, LV doesn't like loose ends
 - so pick up the open-ended wire and connect it to the sane one: bugger! now everything is messed up... Numeric
 - LV comprehensively checks while you code and lets you immediately know if something is wrong



1.231

1.23




Introduction to LabVIEW

Stephan Eisenhardt

And More Wire Skills

□ Of course, we know what is wrong:

- mouse-over a segment of the offending wire
- this time precisely, for the pointer being an arrow click once to select a segment double-click to select up to the next branches triple-click to select everything up to the I/O ports
- get rid of only the offending part, in either way:
 - Ctrl-x
 - right-click: Delete Wire Branch



- now select one of the remaining segments and drag it around
- try different positions and different segments and observe:
 - how does it look if two wires cross each other? look closely...
 - what happens if you try to produce a loop in one wire?
- note: the actual connection to the elements is in the centre of the I/O ports, i.e. under the icons

Introduction to LabVIEW

Stephan Eisenhardt

Wiring needs Order


- □ Is your wiring well messed up making the code incomprehensible?
 - right-click: Clean Up Wire
 - ...aaah... now I get it




- only: the LV algorithms to order wires are not perfect/ergonomic, at times annoying, and in a few cases do bad choices (this may or may not have improved now since v7.1), expect that:
 - at wiring the layout may differ depending on which direction you go between the ports of two elements
 - 'Clean Up Wire' nearly always will come up with a different layout than the wire placement
 - labels are not cared for: wires run straight through them
 - in some circumstances (usually with bigger structures) the wires may leave a port in the opposite direction than one would expect, and emerge under the other side of the object (especially hard to decipher in third-party code...)
- plea: use manual re-ordering where needed, improve readability of code


Finalise the Sum: the Indication

- You want to display the calculated sum. What is the best way think:
 - I suggest: right-click the output port of the sum and create an Indicator
 - this way a new variable is allocated, already wired and it inherits the necessary properties (e.g. data type, name label of the output port)
 - and where did the data display end up?
 - here is a short-cut: right-click the new element \rightarrow Find Indicator
 - and the way back: right-click the new element \rightarrow Find Terminal
 - note the different features of the Indicator:
 - greyed data display, no controls
 - thin-lined terminal, with data entering from the left
 - though you still can edit everything as for the terminal, try:
 - radix, Unit Label, data value, make some data default (even if it doesn't make sense at this stage...)
 - now change the data representation to U16 and enter (x,y)=(2.5, 3.1)
 - run the VI: File Edit View Project




Introduction to LabVIEW






Selecting

- □ You have messed up the configuration? OK, start again:
 - on the FP click the indicator and remove it (Ctrl-x)
 - notice that the 'Run' button breaks, ngrrr... :
 - go to the BD and look what is wrong: aha, a loose wire was left over
 - right-click the wire stub \rightarrow Remove Loose Ends: OK
 - could we have done that easier? yes!
 - 2x Ctrl-z to revert to the problem...
 - on the BD now point to the background and lasso the indicator and wire..
 - selecting this way needs a bit of practise, as:
 - all elements (or labels) touched will be selected
 - wire segments only will be selected if >50% covered
 - hold Shift, point & click to select and deselect other items







⊻±⊻ı



Programme Execution

File Edit View Project

and

□ Have you noted?


- while you code the 'Run' button changes between
- Elle Edit View Project, i.e. enabling and disabling the execution of the VI
- the LV interpreter is constantly checking whether your VI complies with all rules
- some basic reasons for broken VIs you know already:
 - at least one input port of one element does not get data
 - unused wire stubs are present, i.e. data is not properly delivered (though it is OK to leave output ports unwired, i.e. unused)
 - mismatch in data type between output and input of two connected elements (though the input port of an element may provide limited data coercing with default casting rules – it is not recommendable to rely on this)
- a working 'Run' button means that your program compiles and executes, this does not state anything about the programme logic, i.e. whether it makes sense... that remains your job...

Saving



- it's now a good time to save your VI,

you may have to return to it




- Ctrl-s to select path and name of the VI (or to overwrite if existing)



- the '*' always indicates unsaved changes, even if only of cosmetic nature
- and you have now named the VI
- if you exit a VI with unsaved changes, LV will ask you whether they should be saved
 - if there are any sub-VIs with unsaved changes which this VI depends on, LV will also ask you whether you want to have them saved as well
- use Ctrl-Shift-s if you want to save your current VI and all it depends on
- use menu File \rightarrow Save As for more options

Various Manifestations of Elements


- Elements on the BD and FP can have various manifestations:
 - on the BD: terminals and indicators can be shown as
 - icon
 - variable
 - with exactly the same functionality
 - to toggle: right-click \rightarrow View As Icon




Numeric 0

- on the FP: controls and indicators have a wide range of representations
- to select style: right-click \rightarrow Replace \rightarrow Num Ctrls \rightarrow <Style>
 - some examples of inputs for double variables:
- functionality can be configured to suit the task:
 - change the input range of the slide: double-click the axis labels to edit
 - what happens if you edit a centre axis label?
 - operate the knob as input for a double, now change the representation to U8
 - how does the operation of the knob change?
 - right-click a slide/knob \rightarrow Mapping \rightarrow Logarithmic

– note: you cannot enter "0" anymore in DBL representation, but in U8
 Introduction to LabVIEW
 Stephan Eisenhardt






Debugging: Coding Errors

- Now we deliberately break the VI to learn how to investigate Numeric problems:
 - convert the indicator into a terminal: on the FP or BD right-click the indicator \rightarrow Change to Control
 - explore the ways to find out about the error:
 - in the BD mouse-over the red cross/broken wire to get a quick explanation/suggestion
 - hit the broken run button File Edit View Project • to get an Error List




with a more verbose explanation for each error/warning Show Error will focus back to the location on the BD Help will open LV Help with an exhaustive discussion



menu Help \rightarrow Explain Error...

still after resolving problems like the above there may be run-time errors, they usually come with error codes, here you can get them better explained

Stephan Eisenhardt



1.23

x+y





Debugging: Programming Logic

- When LV finds all rules obeyed the run button is unbroken, but this does not check for flaws in the programming logic:
 - these you can sift for in the debugging running mode:
 toggle the running mode with the lamp button
 single or continuous operation then will execute in slow-motion
 - you can enter arbitrary values at the input terminals before execution
 - on the BD you can watch the data packets moving between the elements
 - the values of the data packets are shown at the I/O ports
 - pause III or stop I operation at any point
 - when paused, the element just being executed will blink
 - mouse-over the wires and I/O ports to inspect the values
 - step-buttons allow to execute the code element by element:
 - branch into sub-VI and stop at its first element
 - step over this element
 - finish this VI and stop after its call in the parent VI

Introduction to LabVIEW

Stephan Eisenhardt






Numeric

Numeric 2

13.78

Debugging: using Probes

- Probes are useful to watch values from various locations:
 - probes are updated as the VI executes
 - right-click a wire \rightarrow Probe
 - this will add the probe to the Probe Watch Window
 - double-clicking a probe will focus back to the BD



🎦 🚼 🗙		2	~	Probe Display
Probe(s)	Value	Last Update	~	Numeric
Sum_DBL.vi				15.5
[1] Numeric	15.500E+0	21/09/2010 12:13:00		15.5
[2] Probe	19.280E+0	21/09/2010 12:13:01		

- right-click a probe to access further options, e.g. copy data
- buttons provide further options to manage your probes:
 - opens a new window for the selected probe use it next to the BD or FP during execution
 - selects all probes in the Probe Watch Window use Shift-click and Ctrl-click to alter the selection
 - 🔀 remove a probe from the list
 - collapse/expand the probe display on the right



Introduction to LabVIEW


Stephan Eisenhardt

Debugging: using Breakpoints

- Program execution will pause at breakpoints:
 - right-click a wire or element → Breakpoint → Set Breakpoint the execution will pause when a data packet passes the point and you can take it step-by-step from there
 - to manage the points:
 - right-click a breakpoint → Breakpoint → Clear Breakpoint to remove the point
 - right-click a breakpoint → Breakpoint → Disable Breakpoint to keep but ignore the point
 - multiple points can be easily handled with the Breakpoint Manager:
 - right-click a breakpoint \rightarrow Breakpoint \rightarrow Breakpoint Manager
 - similar to the Probe Watch Window us can use:
 - to enable a selected point
 - o to disable a selected point
 - to remove a selected point
 - selects all points
 - use Shift-click and Ctrl-click to alter the selection

Introduction to LabVIEW

Stephan Eisenhardt




Numeric

VI Name	Object Name	State 🔺 🥚
Sum_DBL.vi	Wire	
Sum_DBL.vi	Add	Č C
		X
		✓ ‡6
	Close	Help

Debugging: using Conditional Probes

- □ An even more powerful tool is the Conditional Probe:
 - right-click a wire \rightarrow Custom Probe \rightarrow Conditional xxx Probe
 - depending on the data type you can set conditions on which this probe

acts as a breakpoint



1 🕄 🗙		2	Probe Disp	olay	
Probe(s) Sum_DBL.vi [1] Numeric [3] Numeric 2	Value 15.500E+0 3.780E+0	Last Update	🗹 Equ	ater than	

 you can switch between the data view and the conditions view in the Probe Watch Window as well as in the new window to permanently

watch the probe:

🔁 [3] Nu	ımeric 2	_ 🗆 🛛		🔁 [3] N	umeric 2		_ 🗆 🛛
Data	Condition			Data	Condition	1	
				Pause if a	ny of the fo	llowi	ing:
				🛃 Equal	to	ч	3.5
a 3.78				🗹 Great	er than	а	3.7
				🔽 Less t	han	а	3.3
			L				

otherwise it behaves like ordinary probes

Keeping Order

- If VIs become more complex keeping order is essential to maintain readability (for others and yourself...):
 - even though LV provides some support for that Image and Image a
 - lasso/select some elements on the DB/FP and play with the tools to see:


. ▼	<u> </u>	\$~ 🛰	1
Horizo	ontal Ce	nters	
[•=	• <mark></mark>		
	• • •		

- Image: Align Objects: align object boundaries or centres of selected elements
 - object boundaries are always the square around the full visible items, i.e.
 depending on the manifestation and including labels, but only if they are shown...
 - when aligning to top/bottom/left/right always the xxx-most will be the reference
 - when aligning to centres alignment will be done to the centre of the selection
- The Distribute Objects: arrange objects apart
 - object boundaries as for alignment
 - distribution takes place within boundaries of selected objects
 - compression makes objects touch
- Main Clean Up Diagram: fully automatic rearrangement
- 🔄 Reorder:
 - group or ungroup selections
 - select stacking order if elements overlap (I advise strongly to avoid overlaps!)

Introduction to LabVIEW

Keeping Order: Examples

□ Here is some taste of what the automatic tools can do:




- so: use the tools with care and don't shy away from manual improvement
- and there is always Ctrl-z ...

Introduction to LabVIEW

Extra Labels & Decorations

- □ Be kind to others and yourself add explanatory labels:
 - you will value them when you revisit your code after some time
 - double-click on the background of the BD/FP to start editing a label
 - double-click the label to re-edit
 - select to move or change the font 13pt Application Font
 - right-click the label for further options




- decorations on the FP or BD like frames, bars or arrows may significantly improve the overview:
 - FP: Controls Palette → System → Decorations
 - BD: Functions Palette \rightarrow Programming \rightarrow Structures \rightarrow Decorations
- colouring them is possible via Property Nodes and Colour Boxes (but that is too complex to be convenient for a 'quick tint' of some elements)

The Interface

□ Encapsulation is mediated via the VI interface:

- lets customise the interface of our example (currently default):
 - double-click icon to open the icon editor:
 - this 32x32 pixels image is helpful if you design modules which you call later
 - you can play with it if you have time, e.g. you can alter it to:
- more importantly, manage the hidden connector layer:
 - on the FP right-click the icon \rightarrow Show Connector
 - in the default layout you see 6 input ports on the left and 6 output ports on the right
 - right-click connector \rightarrow Pattern, to select the layout suitable for this VI
 - or right-click connector \rightarrow Add/Remove Terminal, to adjust left or right side
 - click a terminal in the connector to select note the wire tool appearing and next select an input or output element on the FP to make a connection
 - break such a connection by right-click connector \rightarrow Disconnect This Terminal
 - note the colour of connected terminals representing the data type
 - note that (input: left, output: right) is only a convention, you can cross-connect
 - configure the connection with right-click connector \rightarrow This Connection Is \rightarrow XXX
 - mouse over the connector or icon: note the automatic documentation in the
 - Context Help now your VI is a new building block, ready to be called

Stephan Eisenhardt



_ 0



x+y

LabVIEW



1111

Context Help

Numeric

Numeric

Sum_DBL2.vi

E

_ 0










Configuring VI Behaviour

- To individually configure the behaviour of your VI during edit and running:
 - Ctrl-I or menu \rightarrow File \rightarrow VI Properties, explore categories, e.g.:
 - General: Current revision counter increments every time you save the VI
 - Editor Options: you can alter the grid sizes for FP and BD
 - the defaults are sensible...
 - Window Appearance: customise features available at run-time
 - look up the options in Custom \rightarrow Customize
 - Execution: Allow debugging: switch off if performance matters
 - Execution: Reentrant execution:
 - several instances of VI can run in parallel(second can start before first has finished)
 - use 'Preallocate clone for each instance' if they run independently
 - use 'Share clones between instances' for recursive operation

configure the example this way for later use in a recursive setting (OK), and save the VI



• Execution: Run when opened: e.g. when the user shall not edit...

Introduction to LabVIEW

Configuring LabVIEW Behaviour

□ To configure the general behaviour of LV, across the VIs:

- menu: Tools \rightarrow Options, explore categories, e.g.:
 - Front Panel: Front Panel Grid: Enable Grid Alignment
 - toggle the snapping to the grid if objects are moved
 - Block Diagram: Block Panel Grid:
 - Enable Grid Alignment
 - useful for structured coding
 - Environment: Maximum undo steps per VI:
 - set to 99 (max value, default = 30)
 - Environment: note the default colour scheme

Options		_ 🗆 🗙
Category A	Environment	
Front Panel Block Diagram Controls/Functions Palettes Environment	*Changes to marked options will take effect the next time you start LabVIEW.	
Paths Printing Source Control Menu Shortcuts Revision History Shared Variable Engine V Sarver Web Server Web Services: Security	Colors Use default colors Front Panel Block Diagram Coercion Dets Scrollbar Blink Foreground Blink Background Blink Background Menu Text Menu Background Custom Colors	
	OK Cancel	Help

- Menu Shortcuts: manage you fancy key sequences
- VI Server: to configure conditions under which a remote application can call a local VI

Introduction to LabVIEW

Local Variables

- □ Traditionally LV knows Local and Global variables:
 - the scope of a Local Variable is the VI it resides in:
 - it can be bound ('linked') to any I/O terminal of the VI
 - writing/reading to the Local Variable is as if writing/reading the terminal itself
 - this way one also can write to a control or read from an indicator
 - just to play with it:
 - to add one, right-click BD \rightarrow Programming \rightarrow Structures \rightarrow Local Variable
 - to link it: right-click \rightarrow Select Item \rightarrow xxx
 - change I/O direction: right-click \rightarrow Change to Read
 - and for now get rid of it again...







Numeric

Global Variables

□ Global variables can be powerful and dangerous:

- the scope of a Global Variable is, surprise, global:
 - it resides in a special VI, that only consists of a FP to contain global variables
 - any I/O terminal in any VI can be linked to it
 - several distinct VIs with global variables can be maintained
 - matching names of two global variables from different VIs are possible...
 ... and bear a high risk of confusion (hard to debug!!)
 - they are useful for general configurations
 - they may be used for communication between independently running VIs
- just to play with it:
 - to add one, right-click BD \rightarrow Programming \rightarrow Structures \rightarrow Global Variable
 - to create it: double-click \rightarrow "Global X Front Panel" opens \rightarrow place a Numeric
 - to link it: right-click \rightarrow Select Item \rightarrow xxx
 - note: only the symbol distinguishes the global variable now from the terminal advice: use clear, distinctive naming conventions!
 - to change the Global VI to select from: right-click → Replace → All Palettes
 → Select a VI

Introduction to LabVIEW

Stephan Eisenhardt

Oumeric

Nun	neric
0	

3?

Shared Variables

□ Shared Variables are a relatively new concept in LV:

- they are designed for fail-safe inter-VI communications:
 - e.g. they can be buffered and protected against multiple write accesses
 - but they are only available on Windows and Real-Time platforms
 - inter-communications are supported for:
 - different parts of a BD where wiring is difficult or impossible
 - between independently running VIs on one host
 - between independently running VIs on different hosts on the same network
 e.g. between a master VI on a control PC and an autonomously running VI on a remote Real-Time data acquisition system
- to get a flavour of Shared Variables you may feel adventurous and try:
 - embed the VIs into a project (menu: Project \rightarrow New \rightarrow Add VIs)
 - in the Project Explorer window: select an item and use New to create a shared variable, then configure to need (to complex to discuss here...)
 - ... we don't need all that now: to get rid of all quit LV without saving...

A Recursive Example I

□ We start from our simple example again and build a recursive VI:


- we need to rebuild:
 - open your VI and save it e.g. as Modulo.vi, update the label text
 - rename the controls to 'sum' and 'divisor'
 - remove the 'x+y' indicator
 - add: DBL control 'dividend', U32 control 'iteration'
 - add: DBL indicator 'sum out', U32 indicator 'iteration out'
 - ensure for all controls/indicators that '0' is the default value
 - edit the VI connector to 4 input and 2 output ports and wire up the terminals
 - edit the VI icon to bear the name 'modulo'
 - add a 'x>y' comparison & link the output of the '+' to x and the 'dividend' to y
 - add a Case Structure and link the '?' to the output of the '>' operation
 - on the BD right-click the 'divisor' control, create a local variable and place it in the 'false' case and change it to Read; repeat for the 'dividend'
 - add a '+1' increment in the false case
 - right-click $BD \rightarrow Select VI... \rightarrow place Modulo.vi in false case$
 - wire up the VI inputs/outputs: the sum from the '+', the iteration via the '+1'

Introduction to LabVIEW

A Recursive Example II

□ How does your VI look like? Maybe like this:

- note the 4 tunnels (= feed-through terminals):
 - to hand over data: filled on the left \rightarrow OK
 - broken on the right: because there ins no input in the 'true' case yet...
- finish the VI:
 - in the 'true' case: feed through the 'sum' and 'iteration'
 - multiply the 'iteration' with the 'divisor'
 - from the product create an indicator 'quotient'
 - move the 'quotient' outside the case structure, and wire back to the product
 - in the 'true' case: subtract the 'quotient' from the 'dividend'
 - from the difference create an indicator 'remainder'
 - move the 'remainder' outside the case structure, and wire back ...
 - upgrade the connector to 4 outputs and connect the 'quotient' & 'remainder'
 - fix the link to Modulo.vi: right-click \rightarrow Relink To SubVI
 - fix the 'false' case by wiring the 'quotient' & 'remainder'
- try out the VI with values>0 for dividend & divisor: it should work...



A Recursive Example III

- □ The basic functionality is OK now:
 - here is how my Modulo.vi looks like:


			ools <u>W</u> indo	w <u>H</u> elp						F
	口雪	13pt Applie	ation Font	- 1	· •0= •	₩ ~	⇔ -		3 -	
		sum					sum out			
		(j) 0					12			
		divisor	simple e	ducational for recursion	2			remainder		
		÷)3			2			1		
		dividend						quotient		١.
		÷) 10						9		
		iteration					iteration	nout		
		÷)0					3			
			_		_					
			۲۲	ont	Pa	an	el			
()									>	


– but I am not happy:

there are four major flaws

- 1) the FP shows to much information
- 2) the default values for dividend and divisor (0) cause an infinite recursion
- 3) the VI is not safe against ill chosen input values

4) the VI immediately starts execution: a user dialog for input is desirable Introduction to LabVIEW Stephan Eisenhardt





50

A Recursive Example IV

- □ I also get a bit ambitious and want the VI to have a pop-up GUI:
 - solving issue 1) is easy:
 - for the sum/iteration control/indicators \rightarrow Advanced \rightarrow Hide control/indicator
 - issues 2) & 3) actually can be solved in one go: with a test against
 - but the GUI needs infrastructure and it is sensible to write a wrapper VI:
 - create a Modulo_GUI.vi menu: File \rightarrow New VI (Ctrl-n)
 - add a 'Prompt User' dialog box to the BD:
 - right-click BD: Programming \rightarrow Dialog & User interface \rightarrow Prompt User
 - configure the dialog box to ask for dividend & divisor
 - wire collected data to the tests
 - the output of the tests should 'enable' two Display Messages informing about the numerical requirements for the inputs
 - place a While Loop around everything:
 - change the Loop Condition to Continue If True
 - add an OR ≫ between the two test results and the Loop Condition
 - right-click loop edge to add two Shift Registers a ,wire user values to them
 - viola: you have an input dialog and safe-guard against unsupported inputs

Express VIs

- The User Prompt and Display Message structure you just have used are Express VIs:
 - their main purpose is to make the life easy for the programmer
 - they provide quickly configurable solutions for standard tasks
 - note: corresponding fundamental functions may offer more functionality and control than Express VIs, but are significantly less convenient to use
 - Express VIs are most powerful in the I/O towards:
 - user
 - files
 - (virtual or real) instruments (Signal Express)


A Recursive Example V

□ Continue towards the output GUI:

- first we do it 'per pedes', i.e. with fundamental functions:
 - place the Modulo.vi on the BD
 - wire the user variables and an initial '0' to the 'sum' terminal of Modulo.vi
 - convert the 4 DBL variables to strings:
 Programming → Sting → String/Number Conversion → Number To Fractional String
 - add text constants (Divident:)and Carriage Return Constants
 and use the expandable Concatenate Strings function
 to build a suitable output text
 - add another Display Massage to present the calculated result to the user
 - · for the wiring to look reasonable: use an iterative combination of
 - Clean Up Wire
 - manual corrections

A Recursive Example VI

□ How does your VI look like? Here is mine...




 now replace the string configuration with a Build Text Express VI... Introduction to LabVIEW
 Stephan Eisenhardt


Hierarchies

□ We have a little hierarchy of VIs now: let's have a look at that...

- Modulo_GUI.vi: menu: View \rightarrow VI Hierarchy
- note the difference shown if done in Modulo.vi
- explore the possibilities, e.g.:
 - switch on/off branches

 - include/exclude VI Lib / Global Variables <a>[m]
 - menu: Tools \rightarrow Find VIs on Disk...
 - right-click VIs:
 - Show/Hide All SubVIs
 - Show All Callers
 - Find All instances
 - double-click: to open VIs





Introduction to LabVIEW

VIs as Executables

- Finally: VIs ready for use can be compiled to stand-alone executables:
 - this has benefits in the performance
 - the code ins protected against (accidental) alteration
 - menu: Tools \rightarrow Build Application (exe) from VI...
 - from: Modulo_GUI.vi
 - this generates a Project:
 - Modulo.lvproj
 - Modulo.aliases
 - and opens the VI Properties dialog:
 - select target name and location
 - Build

– it runs...

... as a slim

... without LV

• bingo:

it runs	🗷 Build status 📃 🗆 🔀
without LV open	Modulo_GUI
as a <mark>slim</mark>	The build is complete. You can locate the build at C:\eisenhar\work\UoE\SeniorHonours\builds\Modulo_GUI\Modulo_GUI.
executable	Warnings
	Done Cancel

Build Application from VI	×
Use this dialog box to create a stand-alone application from Modulo_GUI.vi. When you click the Continue button, LabVIEW creates a new build specification in a new project and diplays the Application Properties dialog box. To build the application, click the Build button in the properties dialog box. LabVIEW project (.lvproj) to create	
C:\eisenhar\work\UoE\SeniorHonours\LabVIEW_introduction\Modulo_GUI.lvproj	
Continue Cancel Help	j

Modulo_GUI2 Properties		
Category Information	Information	
Source Files Destinations Source File Settings Icon	Build specification name	
	Modulo_GUI	
Advanced Additional Exclusions Version Information	Target filename Modulo_GUI.exe	
Windows Security Shared Variable Deployment	Destination directory	
Run-Time Languages Preview	C:\eisenhar\work\UoE\SeniorHonours\builds\Modulo_GUT\Modulo_GUT	
	Build specification description	
	Build OK Cancel	Help

Introduction to LabVIEW

No Coding Experience Yet?

- □ You are still worried that you have no coding experience yet?
 - you have!
 - just the set of functionality is limited yet...
 - but now you know a comprehensive set of tools and practices to find you own way around LabVIEW
 - it is time to make use of them
 - the planes of LabVIEW are vast and fruitful...
 - to sustain yourself you have been given bow, arrows and a knife
 - and you learnt how to use them
 - now go: hunt & explore

Best Coding Practices


- □ But it's a long way to write good, complex code...:
 - ... the following is not needed for small, quick hacks but becomes beneficial for bigger projects
 - in LV Help: search for 'best practices', e.g. look at:
 - LabVIEW Style Checklist: revisit as you progress gaining experience

Creating Documentation

>

- VI Memory Usage: dive deep...
 - get aware of optimisations
- VI Execution Speed:
 - performance-orientated...
- in LV Help: use Locate
 - e.g. look at:
 - \rightarrow Fundamentals
 - \rightarrow Development Guidelines
 - \rightarrow Concepts
 - \rightarrow Common Development Pitfalls

Introduction to LabVIEW



difficult it becomes to avoid overlap.

who can work on a project effectively. With fewer people, there is less overlap. You can <u>partition</u> the project so each person works on a particular section. The more developers you add, the more

Recommended Resources I

- □ Helpful resources from NI are:
 - LV menu: Help \rightarrow Web Resources... (the official path)
 - \rightarrow Technical Resources \rightarrow Getting Started with LabVIEW: then e.g.
 - → Learn LabVIEW Basics
 - \rightarrow On-Demand LabView Training (registration/login required)
 - NI: LabVIEW Introduction Course Six Hours (the traditional approach)
 - http://zone.ni.com/devzone/cda/tut/p/id/5241
 - NI: Support Forum, Technical Support & Developer Zone (the Full Monty...)
 - http://forums.ni.com/ni/
 - http://www.ni.com/support/
 - http://zone.ni.com/dzhp/app/main
 - NI: LabVIEW Object-Oriented Programming (available only from LV 8.5)
 - http://zone.ni.com/devzone/cda/tut/p/id/3574
 - LabVIEW mailing list: read by experienced/helpful people, some from NI
 - http://www.info-labview.org/

Introduction to LabVIEW

Recommended Resources II

- □ Selected unofficial web and offline resources are:
 - LabVIEW Wiki: with LabVIEW Turorial
 - http://labviewwiki.org/Home & http://labviewwiki.org/LabVIEW_tutorial
 - Rensselaer Polytechnic Institute: Tutorials in G (for first practical steps)
 - http://www.cipce.rpi.edu/programs/remote_experiment/labview/
 - Connexions: online course on LabVIEW Graphical Programming (LV 7.1)
 - http://cnx.org/content/col10241/1.4/
 - LAVA: LabVIEW Advanced Virtual Architects (for advanced problems)
 - http://lavag.org/
 - Univ. of Utah, Dep. of Physics & Astro.: LabVIEW teaching materials
 - http://www.he-astro.physics.utah.edu/~jui/3620-6620/
 - \rightarrow Supplement Materials
 - \rightarrow Lectures and accompanying Materials
 - Stackoverflow: Q&A forum for programmers \rightarrow search for LabVIEW
 - http://stackoverflow.com/questions/tagged/labview
 - Books: LabVIEW for Everyone & The LabVIEW Style Guide

Stage 2 : More Useful Functionality

Senior Honours Projects

- □ You have now a general overview of the functionalities of LV:
 - not all of the discussed you will need to complete your Senior Honours project
 - and there is more functionality you should familiarise yourself with before you start designing and coding for the project



– T:



– T:

Conclusions

